

Markers

2019 PRODUCT GUIDE

REFERENCE

"Each radiograph must include an appropriate marker that clearly identifies the patient's right (R) or Left (L) side. Medicolegal requirements mandate that these markers be present. R and L markers must be used with CR and DR digital imaging. The development of digital imaging and the use of CR and DR have enabled an environment in which the R and L markers can be placed on the image electronically. This is not recommended because of the great potential for error and legal implications. Anatomic markers should be placed on the CR cassette or DR table similar to screen-film cassettes." [Merrill's Atlas of Radiographic Positioning and Procedures, 9780323073219, text extract]

Markers can be made with letters or numbers or a combination of both.

SOLUTIONS FOR EVERY IMAGING NEED!











Initials Numbers

2 Lines of Letters

2 Lines of letters/numbers

NFW Marker Mates

A new, reusable sticky strip that adheres markers to receptors or equipment. Simply peel off the protective film and apply to the back of your markers. Each pair of Marker Mates will last approximately 30-days under normal use and care. No more need for tape when you can instantly mate your markers!

NEW AND IMPROVED

- The world's strongest x-ray adhesive
- 3x Stronger than standard marker adhesives
- Tested up to 6 months of continuous use in clinical studies
- Works with every type of marker!
- Easy to wash, clean and resue
- Reduces the spread of infections and disease

#MM-XRAY #MM-XRAY10 #MM-XRAY25 1 pair 10 pair 25 pair

Save money and order in bulk today!

#MM-XRAY100 #MM-XRAY250 #MM-XRAY500 100 pair 250 pair 500 pair



Elite Style

Elite markers have the most DURABLE aluminum casing. These markers provide the absolute best in QCI. Easily identifiable. RED for "R" markers and misc. words / BLUE for "L" markers.





#1E ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#4E set



#7E ea. set
R & L - 1"
2 initials - 3%" or 3 initials - 1/4"
Without initials, order:
#5E set





#11EN ea. set

Neonatal Markers R & L is 3/16" and comes with 2 initials or numbers that are 1/8".

The overall size of the markers is 5/8" tall and 3/4" wide.





#8E ea. set
NEW Elite PLUS
R & L - 5/8" Comes with
2-3 initials per line - 1/4"






#TA-137 ea. set R & L is 1/2". 2 initials - 1/4" NEW Without initials, order: #TA-18 set





#2E ea. set
RIGHT & LEFT - 1/2". 2-3 initials - 1/4"
Without initials, order:
#122E set



#6EB ea. set
R & L - 1/2", arrow - 1". 2-3 initials - 1/4"
Without initials, order:
#6EA set



#A-3 each Arrow - 1" Elite Style

Hundreds of custom combinations available!

Identifier Style

Identifier markers have high-impact polystyrene casing. Enhances quality control with fewer repeat exams. Provides basic, yet personalized image identification.





#TA-1 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-4 set





#TA-7 ea. set
R & L - 1"
2 initials - 38" or 3 initials - 1/4"
Without initials, order:
#TA-5 set





NEW #TA-8 ea. set R & L is 5/8", and comes with 2 or 3 initials per line- 1/4".





#TA-2 ea. setRight & Left - 1/2", 2 or 3 initials - 1/4"
Without initials, order: #TA-122 set



#TA-LK ea. set Look Right / Look Left Arrow + "LOOK" - 1/2"





#TA-9 ea. set R & L is 1", and comes with up to 3 initials 3/16.





#TA-002 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-001 set



#TA-A3 each Arrow - 1" Identifier Style

Hundreds of custom combinations available!

Position Indicator Markers

Beads Indicate Patient Position

Our MOST popular Position Indicator Marker available in a variety of colors. Mix and match to make your work day fun! \$32.00 ea. set with initials

















#TA-13 ea. set Elite Style R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4" Without initials, order: #TA-14 set

#TA-13-CG

#TA-13-PK

#TA-13-BK













#TA-13-GPK

#TA-13-TR

#TA-13-CY

#TA-13-CO

#TA-13-PW







NEW #TA-11 ea. set Elite Style R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4" Without initials, order: #TA-12 set



#TA-17 ea. set Identifier Style R & L - 7/8" Position 2 initials - 3/8" or 3 initials - 1/4"





#TA-16 ea. set Identifier Style R & L - 7/8" Position No initials

Personal Touch Colors & Shapes

Make your mark with our customizable styles!

Sold as a set

R & L - 5/8"

2 initials - 3/8" or 3 initials - 1/4"



















BLACK SOLID #PT-BK

BROWN SOLID #PT-BR

GREEN SOLID #PT-CG

LIGHT BLUE #PT-PW















YELLOW SOLID #PT-CY

ORANGE SOLID #PT-CO

PURPLE SOLID #PT-CP

PINK SOLID #PT-PK

Mix and Match Your Favorite Colors!

Want to add some "fun" to your work day? Mix and match colors. It's easy! Choose the color from above and "L" or "R".

For example: (for technologist Terri Andrews)

PT-BK/L is a Black Left Marker with initials "TA"

PT-CO/R is an Orange Right Marker with initials "TA"





PT-BK/L

PT-CO/R

Personal Touch Glitter

Time to make your mark with our unlimited custom combinations!













RED/BLUE GLITTER #PT-GR ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

PURPLE GLITTER #PT-GP ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

LIGHT BLUE GLITTER **#PT-GB** ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

YELLOW GLITTER #PT-YG ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"











NEW WHITE GLITTER #PT-GWH ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

PINK GLITTER #PT-GA ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

ORANGE GLITTER **#PT-GO** ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"

GREEN GLITTER #PT-GG ea. set R & L - 5/8" 2 initials - 3/8" or 3 initials - 1/4"







SKULLS #PT-SKW ea. set





SKULLS #PT-SKW-BO ea. set

Image Scaling Marker Device

Scaling a digital image is mandatory in order to assure the highest level of diagnostic accuracy. It is critical to know the difference between the size of the anatomy shown on a digital monitor as compared to the actual anatomy. An easy to use, easy to clean, easy to understand calibration tool is a must, and a Rollr Precision Marker is the answer!

The patent pending radiolucent sphere holder and a size indicator sphere come together to allow perfect visibility of the sphere on the image as taken from any direction. A millimeter marker on the sphere holder indicates the size of the sphere that is being used.

The suction cup base allows the unit to be firmly secured to a flat table or a vertical wall bucky board that is close to the patient. The small footprint of the suction cup does not interfere with the patient, the technologist, nor limit its use.

The flexible gooseneck arm is perfect for positioning the calibration sphere as close as possible to the critical anatomy. The arm is covered with a smooth, crevice-free surface that allows for easy sanitation.

The autoclavable unit comes with the suction stand and one calibration sphere. Additional spheres can be bought separately in packs of six (6).



Full unit with 25mm sphere

#A25DR Single

Full unit with 25.4mm (1-inch) sphere

#A24DR Single.

Full unit with 30mm sphere

#A30DR Single




Replacement 25mm sphere (set of 6) #S25DR Set

Replacement 25.4mm sphere (set of 6) #S24DR Set

Replacement 30mm sphere (set of 6) #S30DR Set





CM Ruler Markers

The personalized R & L marker sets or the stand alone 3cm ruler are specifically designed for Teleradiology, CR and DR Imaging, where the image-size does not always match the life-size anatomy. Place a marker on the Receptor Plate and capture the radiographed anatomy. The image of the unique 3cm ruler gives an immediate visual and measureable cue as to the size of the screen image compared to the life size anatomy. Each marker has a lead 3cm ruler. Color coded: L-Blue, R-Red. High density lead characters.


Elite style R&L Markers with patient indicator and 1 centimeter ruler for immediate visual and measurable cue to the size of image captured.





NEW #TA-40DR ea. set Ruler Style R & L - 5/8" 2 initials - 3/8"

NEW #TA-50DR ea. set Ruler Style R & L - 1/2" 3 initials or numbers - 1/8"



#1DR ea. set Ruler Style R & L - 3/4" 2 or 3 initials - 1/4"



#2DR ea. set Ruler Style R & L - 1/2" 2 or 3 initials - 1/4"




#3DR each Ruler only

Wafer-Thin Digital Imaging Markers

Super Thin for SlimBucky Clearance

Perfect for equipment (CR, DR, etc.) with very thin bucky clearance. Exceptionally thin markers: <3/32" (2mm). Copper filtered for a broad range of KV. Elite style and color-coded.




#WTS-DR ea. set Ruler Style R & L with 3cm spacing



Radi-O-Paque Extremity Rulers

Sturdy 11/4" wide x 1/8" thick. High density lead letters. Brilliant YELLOW color. Spaced in 0.5cm increments.



Wafer Thin

Super Thin for SlimBucky Clearance

Perfect for newer equipment (CR, DR, etc.) with thin bucky clearance. Exceptionally thin markers: <3/32" (2mm). Copper filtered for a broad range of KV. Elite style and color-coded.



#WTS-6B ea.set R & L - 36" 2 initials - 3/16" Without initials (with arrow), order: #WTS-6A set



#WTS-6C ea.set
R & L - 38" 2 initials - 3/16"
Without initials, order:

#WTS-6D set

MBM



#WTS-1 ea.set R & L - 3/8" 2 initials - 3/16" Without initials, order: #WTS-4 set



#WTW-A3 (1/2" arrow)
Each

Oncology, CT and CR/DR

Super Thick for High Radiation

1/4" super thick, high-density lead letters. Designed for Oncology and Industry, and other applications. Elite style markers, color-coded and Techno-Aide tough. Vinyl and expoxy encased for extended life!

WARNING:

The Oncology markers are 1/4" thick and will NOT fit in most low clearance bucky trays.



#R0-107B ea.set R & L - 3/4" 2 initials - 3/8" Without initials, order: #R0-107A set

cat.# (3/4")	description
ROW-01	single letter
ROW-23	2-3 letters
ROW-45	4-5 letters
ROW-67	6-7 letters
ROW-80	8+ letters



Ex.: ROW-23 3/4" letters (1/4"thick)

WORD CHOICE EXAMPLES

2-3 Let	tters	4-5 Lett	ers	6-7 Letters	;	8 or more Lette	ers
AP	LP0	FLAT	ERECT	AP ERECT	SUPINE	ANTERIOR	POST EVAC
PA	MAG	FLEX	HOURS	CLOSED	DELAYED	COMPARISON	POST VOID
ANT	MIN	INSP	PRE OP	INT ROT	FLEXION	DECUBITUS	SEMI ERECT
EXP	OBL	LEFT	PRONE	MOBILE	LATERAL	EXTENSION	STANDING
EXT	RAC	OPEN	PUREE	P EVAC	MINUTES	INFUSION	VALSALVAR
INT	LOA	POST	RIGHT	POST OP	OBLIQUE	INJECTION	VERTICLE W/WEIGHT
KUB	25	TOMO	SCOUT	RECUMB	SUNRISE	INVERSION	WEIGHT BEARING
LA0	600	DECUB	60 MIN	STRESS	UPRIGHT	PORTABLE	WT BEARING

Copper Filtered Filtered for High KV Exams

Our copper filtered markers are specifically engineered for High KV. Crafted with high density lead letters. Available in ELITE or IDENTIFIER. Color-coded, impact resistant.



#COP-1 set ELITE (2 sets of 1E markers: 1 set with copper filter, 1 set without.)



#COP-2 set ELITE R & L - 1/2", 2-3 initials - 1/4" (Filtered & unfiltered on the same marker.)



#TA-CO2 set R & L - 1/2", initials 1/4" - 1/2" Marker includes copper filtration



#TA-C01 set (2 sets of TA-1 markers: 1 set with copper filter, 1 set without.)

Copper Filter Add-On Available! Copper filtering add-on is an additional charge per set

Copper filtering can be added to any marker type for an additional fee. Copper filtering is specifically engineered for High KV scans. Mix and match colors and types, then simply add copper filtering to your customization. cat.#Copper



Standard Thin Clipper

Constructed of high strength Kydex™. Extra thin profile, double-sided,11/8" wide. R and L: 3/4", with or without initials. 2 Initials: 3/8"; 3 initials or words: 1/4".

#CPP-40 (with initials) **#CPP-42** (no initials)





Elite Clippers

Both R & L are on one double clip. Durable Stainless Steel clips. Color-coded: R side-Red / L side-Blue.

Double Personalized

#CP-35 each (2-3 initials)





Portable & Position Double Clippers

#CP-13

(with initials)

#CP-14 (no initials)



#CP-36 (with initials)

#CP-37 (no initials)



Gridcap Thin Clipper

Designed to fit DR, CR receptors with gridcap. Constructed of high strength Kydex TM . Extra thin profile, double-sided,178" wide. R and L: 3/4", with or without



Stencil Clipper

High quality plated steel construction. Stenciled, non-personalized right and left identifiers



Word Choices, Pricing and Accessories

Wafer Thin

(1/4" letters x <3/32"thin)



cat.#	description
WTW-01	single letter
WTW-23	2-3 letters
WTW-45	4-5 letters
WTW-67	6-7 letters
WTW-80	8+ letters

Oncology

(3/4"letters x 1/4"thick)



cat.#	description
ROW-01	single letter
ROW-23	2-3 letters
ROW-45	4-5 letters
ROW-67	6-7 letters
R0W-80	8+ letters

Wall Boards

Walnut grain w/brass hooks.



cat.#	description
HK-32	32 Hook Board
HK-16	16 Hook Board
HK-8	8 Hook Board
HK-50	50 Hook Board

ELITE Words (aluminum casing)





description	example	cat.#	price	cat.#
single letter	3	2E-Y	\$9.00	2E-Z
2-3 letters	EXT	3E-Y	\$11.00	3E-Z
4-5 letters	25 MIN	4E-Y	\$13.00	4E-Z
6-7 letters	STRESS	5E-Y	\$15.00	5E-Z
8+ letters	PORTABLE	6E-Y	\$18.00	6E-Z

IDENTIFIER Words (polystyrene casing)



example 3

EXT

25 MIN

STRESS

PORTABLE

description

single letter 2-3 letters

4-5 letters

6-7 letters

8-15/8-10



cat.#	price	cat.#
2TA-Y	\$8.00	2TA-Z
3TA-Y	\$10.00	3TA-Z
4TA-Y	\$12.00	4TA-Z
5TA-Y	\$14.00	5TA-Z
6TA-Y	\$17.00	6TA-Z

General Marker Set

Elite #GME-60

Identifier #GMI-65



Identifier Set shown (Board: 12" x 15") 3/8" letters / 1/2" numbers



Word Choice Examples

2-3 Let	ters	4-5 Lette	ers	6-7 Letters		8 or more Lette	rs
AP	LP0	FLAT	ERECT	AP ERECT	SUPINE	ANTERIOR	POST EVAC
PA	MAG	FLEX	HOURS	CLOSED	DELAYED	COMPARISON	POST VOID
ANT	MIN	INSP	PRE OP	INT ROT	FLEXION	DECUBITUS	SEMI ERECT
EXP	OBL	LEFT	PRONE	MOBILE	LATERAL	EXTENSION	STANDING
EXT	RAC	OPEN	PUREE	P EVAC	MINUTES	INFUSION	VALSALVAR
INT	LOA	POST	RIGHT	POST OP	OBLIQUE	INJECTION	VERTICLE W/WEIGHT
KUB	25	TOMO	SCOUT	RECUMB	SUNRISE	INVERSION	WEIGHT BEARING
LA0	600	DECUB	60 MIN	STRESS	UPRIGHT	PORTABLE	WT BEARING

Cover Tab Stickers

Ideal for coding mismarked or unmarked film. Color coded: L-Blue / R-Red. Self-adhesive backing. 11/4" diameter with square for initials. Includes 250 stikcers per box. Each box is an idividual letter.

Cat.#	Box	Per Box
TAB-L	1 Box	250 per box
TAB-R	1 Box	250 per box





Marking Pens

Highlights points of interest on a film. Wax markers resist smudging, but are not permanent. NOT to be used in place of lead markers. NOT a legally acceptable marking method.

Order each item per number of packs.

Example: 2ea. #MPY-11

(will be shipped with 10 yellow pens.)

PENS (5 per pack)

cat.#	color
MPY-11	Yellow
MPR-12	Red
MPW-13	White
MPA-19	2-Y, 2-R, 1-W

REFILLS (72 per pack)

cat.#	color
MPY-21	Yellow
MPR-22	Red
MPW-23	White

Refills come in pkgs of 72 pieces.




Skin Markers for General X-Ray, Mammography, CT and Oncology

GENERAL IMAGING: used to identify areas of pain or trauma, soft tissue masses, to document points of entry and exit sites, or mark prominent nipples for identification of nipple shadows on chest x-rays.







Skin Dots

Artifact/nipple identifier. Helps eliminate repeat exams. Lead-free. Easy-Peel tab.




cat.#	BB size
SDM-BB15	1.5mm
SDM-BB20	2.0mm
SDM-BB25	2.5mm
1-5 Boxes	

100 per box

Skin Wings


Artifact/nipple identifier. Extra support for larger, more prominent nipples. Lead-free. Easy-Peel tab.



cat.#	BB size
SDM-W15	1.5mm
SDM-W20	2.0mm
SDM-W25	2.5mm
1-5 Boxes	
100 per box	

Skin Dots (Rad/Onc)

Lead-free artifact identifier. Larger BBs for Radiology and Oncology. Easy-Peel tab.



cat.#	BB size
SDM-BB30	3.0mm
SDM-BB40	4.0mm
SDM-BB50	5.0mm
1-5 Boxes	
50 per box	

Skin Arrow (Rad/Onc)


Lead-free arrow pointer. Highlights areas of interest for Radiology and Oncology. Easy-Peel tab.



cat.#	BB size
SDM-A1	7.0mm
1-5 Boxes 100 per box	

Skin Line/Scar Marker

Lead-free scar identifier. Line contours to complex scar tissue. Flexible adhesive tabs; no scissors required.



cat.#	RR SIZE
SDW-SS612	1.0mm
1-5 Boxes 100 per box (120"/305cm pe	er box)







Economy Nipple-Artifact Markers

Encapsulated 2mm lead ball. Used to mark nipple, lesion, scar, artifact, etc. Neutral color, peel-off, self-adhesive backing. Non-sterile. 100 markers mounted on ten convenient strips. (This product contains lead)



2D/3D Breast Tomo Markers

Used for Digital Breast Tomosynthesis. Low density markers specifically designed to show up clearly when "in slice" and minimize the "slinky" artifact when they're "out of slice."



Mole Marker

Lead-free mole marker. Easy-Peel tab.



box of 100 Hole Size: 1/2"



Palpable Mass Marker

Lead-free Mass skin marker. Easy-Peel tab. cat.# SDM-TSD3

box of 100



Non-Palpable Mass Marker

Lead-free arrow skin marker. Easy-Peel tab. cat.# SDM-ASD3

box of 100



Skin Line/Scar Marker

Lead-free line skin marker. Repeated adhesive tabs allow accurate bending to scar outline. cat.# SDW-SSD3

roll of 120" Line: 1mm

Digital Mammography Skin Markers





Digital Dot

Lead-free nipple identifier. Highlights areas of interest for digital and analog imaging. Easy-Peel tab.

Cat.#

SDM-BB23

1-5 Boxes (100 per box)

CT Skin Markers

BB size 2.3mm



CT Dot: 2.3mm

Lead-free 2.3mm marker. Highlights areas of interest for Radiology and CT. Easy-Peel tab.

Cat.#

SDM-CT1

1-5 Boxes (50 per box)

BB size



CT Dot: 4.0mm

Lead-free 4.0mm marker. Highlights areas of interest for Radiology and CT. Easy-Peel tab.

Cat.#

SDM-CT2

1-5 Boxes (50 per box)



SDG-MSD3

CT Skin Markers

Opaque line skin markers are perfect for CT field edge markings and nearly all CT marking applications. RT Simulation, Treatment Planning and Diagnostic Radiography. Flat marker is perforated every 2" yieldiing 125 markers per box. Lead free and disposable.



cat.# SDW-CT05 (0.5mm) 138 per box



cat.# SDW-CT10 (1mm) 132 per box



cat.# SDW-CT15 (1.5mm) 125 per box



cat.# SDM-RT03 (2mm) 130 per box

Multi-Modality Skin Markers

Provides and accurate reference point for CT, MRI, Therapy, Nuc Med, and PET procedures. Does not indent into the tissue or produce unwanted artifacts. These versatile markers can be used for marking masses, pain location, biopsies, or other areas of interest. Perfect for any modality with a clear accurate reference point. Easy to secure in place. Easy to remove. Manufactured with a medical grade adhesive. MRI safe. Disposable. 50 per box.



CT / MRI

No center hole; appears as bright object on CT and MRI scans. Also highlights areas of interest for radiology and mammography.

cat.# MM3007

1-5 Boxes 50 per box



CT / MRI / THERAPY / RADIOLOGY

2mm center hole suitable for conventional hypodermic needle. Highlights area of interest for Oncology and Radiology.

cat.# MM3005

1-5 Boxes 50 per box



CT / MRI / NUC Med / PET

2mm center hole suitable for injection of radioisotope. Highlights area for Nuc. Med., PET, CT and MRI.

cat.# MM3003

1-5 Boxes 50 per box

MRI Markers

CT / MRI

Easily identifiable. Clearly and accurately communicates the area of interest or concern to the radiologist. Artifact free! Does not indent the skin, preventing a compromised image. 100% Latex-Free and hypoallergenic making it safe for patients. Markers are MR safe.



cat.# SDM-MR2 1-5 Boxes

3/4" capsule 40 per box



cat.# SDM-MR3

1-5 Boxes
1 1/4" capsule

40 per box



Lead Characters

If you have any questions about the size characters you need or are looking for additional sizes and styles, please call our customer service department, or mail us a sample so that we may send you exactly what you're looking for.



Unmounted Lead Characters

style	(3/16")	(1/4")	(3/8")	(1/2")	(5/8")	(1" R&L only)	(1" arrow)
FLAT FACE	XTUF-36 .50	XTUF-14 .50	XTUF-38 .50	XTUF-12 .60	XTUF-58 .60	XTUF-01 1.00	Χ
DEEP BLOCK	Χ	XTUB-14 .90	XTUB-38 1.00	XTUB-12 1.00	(see Flat Face)	Χ	Χ
ARROWS	Χ	XTUA-14 1.00	Χ	XTUA-12 .80	Χ	Χ	XTUA-01 1.00



Flat Face | Deep Block | Arrow

"Loose" lead letters & numbers.

Individual Mounted Characters (Flat Face)

cat.#	lead/plastic size	price	comments
XTSF-3658	3/16" on 5/8"	.80	Fits: LLH-58, LLH-2, LLH-3, Veri-Thin
XTSF-3634	3/16" on 3/4"	.80	Fits: Legal (regular)
XTSF-1458	1/4" on 5/8"	.80	
XTSF-1434	1/4" on 3/4"	.80	
XTSF-1478	1/4" on 7/8"	.80	Fits: Eclipse, Perfection
XTSF-1234	1/2" on 3/4"	.80	
XTIF-3658	3/16" on 5/8"	.80	Fits: Ideal (special offset mounting)
			. 1



Characters are mounted on CLEAR flat plastic tabs to fit a variety of letter holders.

Flat Surface Mounted Words (Flat Face)

description	example	cat.#	price	cat.#	price	cat.#	price	cat.#	price
single letter	3	XVFW-0136	.1.00	XVFW-011	4 .1.00	XVFW-0112	.1.00	XVFW-0101	2.00
2-3 letters	EXT	XVFW-2336	2.00	XVFW-231	4 2.00	XVFW-2312	2.00	XVFW-2301	4.00
4-5 letters	25 MIN	XVFW-4536	4.00	XVFW-451	4 4.00	XVFW-4512	4.00	XVFW-4501	7.00
6-7 letters	STRESS	XVFW-6736	6.00	XVFW-671	4 6.00	XVFW-6712	6.00	XVFW-6701	9.00
8+ letters	PORTABLE	XVFW-8036	7.00	XVFW-801	4 7.00	XVFW-8012	7.00	XVFW-8001	11.00
Lead/Plastic Si	zes	3/16" on 3/4	."	1/4" on 3/4	"	1/2" on 3/4"		1" on 11/2"	



Words made from Mounted Letters glued to the surface of a single CLEAR vinyl strip.

Lead Letter Case

Injection Molded Box. 36 scoop-bottom compartments (for A-Z and 0-9). Side compartment for holder.



Dimensions: 13" x 3" x 13" (1.40lbs)

Ideal Cushion Marker

Double channel with cushions; imprint in the upper section of each channel.

#XSIC-21 each

Holder only. 3 1/2" x 15/16"

Eclipse Cushion Marker

Single channel with cushion; imprint above and below the centerline.

#XSEC-11 each

Holder only. 3 1/2" x 15/16"

Legal Marker

Double aluminum channels: imprint between the two channels.

#XSLG-21 each

Holder only. 3 1/2" x 1 5/8"

Veri-Thin Legal

Double aluminum channels; imprint down and in between the two channels. Great for Bucky Trays!

#XSVL-21 each


Holder only. 3 1/2" x 11 3/16"

Perfection

Single aluminum channel; imprint above and below the channel.

#XSPF-11 each

Holder only. 3 1/2" x 15/16"



Viewable Imprint under Cushion



#XSLG-22 ea set

#XSIC-22 ea set

1ea. holder 2ea. A-Z 4ea. 0-9

2ea. A, E, I, O, U

1ea. 36-letter case (Replacement characters:

Reorder #XTIF-3658)

#XSEC-12 ea set

1ea 36-letter case

Reorder #XTSF-1478)

(Replacement characters: 1/4" on 7/8" plastic tabs.

1ea. holder 2ea. R, L 4ea. 0-9

Complete Set Includes:

3/16" on 5/8" plastic tabs, offset.

Complete Set Includes:

Complete Set Includes:
1ea. holder
2ea. A-Z
4ea. 0-9
2ea. A, E, I, O, U
1ea. 36-letter case
(Replacement characters:
3/16" on 3/4" plastic tabs;
Reorder #XTSF-3634)



#XSVL-22 ea set

Complete Set Includes:
1ea. holder
2ea. A-Z
4ea. 0-9
2ea. A, E, I, O, U
1ea. 36-letter case
(Replacement characters:
3/16" on 5/8" plastic tabs.
Reorder #XTSF-3658)



OLEVELAND CLINIC NAPLES FL

#XSPF-12 ea set

Complete Set Includes: 1ea. holder 2ea. R, L 4ea. 0-9 1ea. 36-letter case (Replacement characters: 1/4" on 7/8" plastic tabs. Reorder #XTSF-1478)

Contact Name	
Address	
City/State/Zip	
Phone	
Email	





	Technologist Name	Initials	Right Color	Left Color	Marker Style
EX.	TERESA QUIRANTE	TJQ	RED	BLUE	TA-1
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					



ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to orders@techno-aide.com

For official Quotes: quotes@techno-aide.com

Any questions can be directed to customerservice@techno-aide.com

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

techno-aide.com