Size Charts

Techno-Aide's **EZ Fit Sizing Tool** quickly and accurately fits people into the perfect Techno-Aide apron size without the need for measuring tapes. However, use these measuring charts to verify your fit if a measuring tape is the only available option.

Techno-Aide standard aprons (front, full wrap, vests, & skirts) are fit using two measurements:

- Width measured at the widest point of the body applicable to apron style
- **Length** measured applicablle to apron style

Visit Techno-Aide.com/ez-fit-sizing-tool for step-by-step measuring instructions and to order your sizing tool today!

FRONT APRONS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Chest Circumference	32-35.5"	35.5-38.5"	38.5-40"	40-46"	46-51"	51-55"	55-60"
Length (from Sternal Notch)	33.5"	35.25"	37.25"	39.75"	41.5"	41.5"	41.75"
Width (at widest point)	20.25"	22.75"	23.5"	24.5"	26.5"	29"	32.5"

EZ FULL WRAP & PARTIAL OVERLAP VESTS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Chest Circumference	32-35.5"	35.5-38.5"	38.5-40"	40-46"	46-51"	51-55"	55-60"
Length (from Sternal Notch)	21"	22"	22.5"	24"	25"	26"	26"
Width (at widest point)	17.5"	18.5"	19.5"	21"	24.5"	27"	29"

EZ CONTOUR & PARTIAL OVERLAP CONTOUR VESTS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Chest Circumference	32-35.5"	35.5-38.5"	38.5-40"	40-46"	46-51"	51-55"	55-60"
Length (from Sternal Notch)	21"	22"	22.5"	24"	25"	26"	26"
Width (at widest point)	17.5"	18.5"	19.5"	21"	24.5"	27.5"	29.5"

EZ REVERSE FULL WRAP VESTS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Chest Circumference	32-35.5"	35.5-38.5"	38.5-40"	40-46"	46-51"	51-55"	55-60"
Length (from Sternal Notch)	18.5"	21.5"	22"	23.25"	23"	24"	23.5"
Width (at widest point)	21"	20"	23"	24.25"	26.5"	29.5"	31.5"

FULL WRAP SKIRTS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Waist Circumference	25-27"	28-30"	31-33"	34-35"	36-40"	41-43"	44-46"
Length (top to bottom)	19"	19.5"	20"	21"	22"	22.5"	22.5"
Hips Circumference	31"	38"	39"	40"	41"	46"	50"

FULL WRAP APRONS

Measurement	Petite	Small	Medium	Large	X-Large	2XL	3XL
Chest Circumference	32-35.5"	35.5-38.5"	38.5-40"	40-46"	46-51"	51-55"	55-60"
Length (from Sternal Notch)	36"	37"	38.5"	40.5"	43"	44"	44"
Width (at widest point)	18.5"	19.5"	20.5"	22.5"	25.5"	28.5"	30.5"



LAP GUARD HALF APRONS

Measurement	Petite	Small	Medium	Large	X-Large	2XL
Width (at widest point)	10"	12"	20.5"	24"	24"	36"
Length (top to bottom)	10"	12"	16"	17"	18.5"	19.5"

BLANKET GUARDS

Measurement	Small	Medium	Large	X-Large
Width (side to side)	18"	24"	24"	24"
Length	18"	24"	36"	48"

PREGNANCY FRONT APRONS

Measurement	Standard	Full	
Width (at chest)	30"	34"	
Length (from shoulders)	43"	45"	

GONADAL GUARDS

Measurement	Small	Medium	Large	
Width (side to side)	4"	6"	7.5"	
Length (top to bottom)	6"	7"	9"	
Strap Length	8-61"	8-61"	8-61"	

THYROID GUARDS

Measurement Standard		Mammo	Pediatric	
Length (strap to strap)	24"	24"	16"	
Bib Height	5.75"	4.75"	4"	

DENTAL DRAPES

Measurement	Child	Adult
Width (at widest point)	18"	23.5"
Length (from neckline)	23.7"	28.4"
Back Flaps (from shoulders)	7.5"	11.5"

PANORAMIC DENTAL DRAPES

Measurement	Adult
Front Width (at widest point)	24"
Front Length (from neckline)	26"
Back Width (at widest point)	25"
Back Length (from neckline)	15.5"

PREGNANCY HALF APRONS

Measurement	Standard
Width (at hips)	27"
Length (from waist)	27"

BREAST GUARDS

Measurement	Child	Adult
Width (side to side)	5"	6"
Length (top to bottom)	6"	7"
Strap Length	8-36"	8-36"

PEDIATRIC FRONT APRONS

Measurement	Standard
Width (at waist)	18"
Length (from shoulders)	25.68"

SCOLIO GUARDS

Measurement	Standard
Width (side to side)	8"
Length (top to bottom)	26"
Neck Strap	8-36"
Waist Strap	49"

SIZING TOOL VIDEO

